

Training Images

These images may be used to assist in training clinicians in the Swift Skin and Wound application.



Step 1.

Using Image #1 the clinician should open a "New Wound" evaluation, capture the image and complete the trace, refine and related documentation.

Step 2.

Using Image #2 the clinician should complete a "Re-Eval" of an existing wound, capture the image and complete the trace, refine and related documentation.



Step 3.

Using Image #3 the clinician should "Resolve" or close the wound. The clinician should complete a "Re-Eval" of an existing wound, after image capture the user should choose "Closed", then proceed to the progress tab to "Resolve" the wound.